



Alcohol Swabs with 2% Chlorhexidine and 70% Alcohol

Disposable

Alcohol swabs with 2% chlorhexidine and 70% alcohol are cleaning wipes that are typically used to sterilize or disinfect surfaces. Chlorhexidine is an antiseptic solution that is effective against a wide range of microorganisms. It is often used in combination with alcohol, which is a strong disinfectant, to create a powerful solution for cleaning and sterilizing. The 2% chlorhexidine and 70% alcohol in these swabs is the concentration of the solutions in the wipes.

Intended Use

It is commonly used for sterilizing surfaces or tools before surgery and other medical procedures, or for cleaning and sanitizing the skin before injections or other medical treatments.

Specifications

- Item: Alcohol swabs
- Material: Non-woven 50g/ m²
- Content: 2% chlorhexidine and 70% alcohol
- Function: Disinfection
- Pack Size: 50*50mm
- Pad Size: 30*65mm
- Feature: Aluminum film bag 73g/ m²

Packing Details

- Package: 100pcs/box,100boxes/ctn
- Carton size: 42.5*28.5x28.5cm
- G.W/N.W: 13/12kg





Instructions

1. Clean the area to be treated: Remove any dirt or debris from the surface or skin that you want to clean and disinfect.
2. Open the packaging: Open the foil pouch or plastic container that contains the swab and remove it.
3. Apply the swab: Gently wipe the area to be treated with the swab. Be sure to cover the entire surface area and pay attention to hard-to-reach areas.
4. Allow the surface to dry: After cleaning and disinfecting, allow the surface or skin to air dry. Avoid touching the area until it is completely dry.
5. Dispose of the swab: After use, dispose of the swab in an appropriate manner, usually in a container designated for biohazard or sharps.

Caution

- External use only: These swabs are intended for external use only and should not be ingested or inserted into body cavities.
- Flammability: The alcohol content of these swabs can make them flammable, so it is important to keep them away from heat sources, sparks or open flames.
- Irritation: The swab should be used with caution on individuals who have sensitive skin or a history of allergic reactions.
- Avoid contact with eyes: If contact with eyes occurs, rinse with copious amount of water and seek medical attention if necessary.
- Follow the instructions: Always read the instructions on the package or follow the instructions given by the healthcare professional for the specific product you are using and the intended use.
- Dispose of the swab properly: After use, dispose of the swab in an appropriate manner, usually in a container designated for biohazard or sharps.
- Not for deep or heavily infected wounds: These swabs are not intended for use in deep or heavily infected wounds or for use on mucous membranes. They should only be used under the guidance of a healthcare professional and following the instructions on the package or given by the healthcare provider.

Storage and Transportation

1. The device shall be stored at a relative humidity of less than 80%, the storage environment shall be kept non-corrosive gases, dry, avoiding sunlight, well-ventilated conditions, and clean.
2. The device shall be taken care during transportation and handling.

