



Bactericidal Wipes

70% isopropyl alcohol

Bactericidal Wipes is produced by non-woven and 70% isopropyl alcohol which has slightly irritating smell. The product is packaged with the Plastic bucket. One wipes saturated with 70% isopropyl alcohol. Single use only which is safe and continent. The product will not cause adverse stimulation to human skin and mucosa, will not cause allergic reactions and other damage, and will not precipitation toxic substances either.

Intended use

Used to disinfect hard surface of non-invasive medical device.

Specifications

- Item: Bactericidal Wipes
- Material: Non-woven 35g/ m²
- Content: 70% Isopropyl Alcohol
- Function: Disinfection
- Size: 14x20cm
- Packing method: Plastic bucket

Packing details

- Package: 100pcs/bucket ,20bucket/ctn
- Carton size: 51x41x19cm
- G.W: 8.5kgs

