



健康



Cryo Vials

Sterile

Cryogenic vials are designed for storing biological material, human or animal cells, at temperatures as low as -196°C but should be used only in the gas phase of liquid nitrogen.

Small, capped vials designed to withstand ultra-low temperatures; various capacities; may be sterile, non-sterile, and/or autoclavable; includes internally and externally threaded vials; ideal for sample storage in vapor state liquid nitrogen.

Specifications

- Item: Cryo Vials
- Further description: External/Internal Threaded Head Cover, Self Standing
- Material: Polypropylene
- Capacity: 1.0ml
- Temperature: $-196^{\circ}\text{C} \sim 121^{\circ}\text{C}$
- Diameter Outer: 12.26-13.79mm
- Length: 40.11mm



Packing details

- Package: 100pcs/bag, 20bags/ctn
- Carton size: 45*32*20cm
- Weight: 5.5kgs

