



Emergency Blanket

Space blanket (known as Mylar blanket, emergency blanket, etc) is an especially low-weight, low-bulk blanket made of heat-reflective, thin, plastic sheeting. The design reduces the heat loss in a person's body, reflects almost 90% of body heat.

Intended use

Used to first-aid devices, survival emergency for warm.

Specifications

- Item: Emergency Blanket
- Material: Aluminum Foil PET
- Content: Blanket,Label,Ziplock Bag
- Function: Survival rescue,keep warm
- Size: 54" x 82"(137CMx208CM)
- Feature:
 - *Reflects almost 90% of body heat
 - *Water and wind-proofed
 - *Remains elastic below freezing
 - *Small enough to fit into your shirt pocket



Packing details

- Package: 1pcs/bag,250pcs/ctn
- Carton size: 44*24x24.5cm
- G.W: 14kgs

