



# Eye Pad


# Disposable

## Applicable for

Operating room, outpatient surgical treatment room and other postoperative wound care and family, personal general wound care.

## Specifications

- Item: Non-woven Eye Pad Dressing
- Material: Spunlace Non-woven : 50g/ m<sup>2</sup>  
Release paper: 80g/ m<sup>2</sup>  
Absorbent pad (cotton) : 150g/ m<sup>2</sup>  
Hot melt adhesive : 32g/ m<sup>2</sup>  
PE film 17g/ m<sup>2</sup>
- Color: White
- Size: 6\*8cm (3.6\*5.6cm)
- Sterile: EO sterile
- Type: Dressings and Care For Materials
- Shape: Oval
- Advantage: High absorbent and softness, easy to use



## Packing Details

- Package: 1pc per sterilization pouch.
- Carton Size: 50pcs/box, 40boxes/ctn  
52\*52\*32cm
- G.W/N.W: 8/7kgs



## Direction for use

- Wash hands;
- Gently clean and dry the eyelid and surrounding area;
- Remove pad from wrapper;
- Apply pad and press down edges firmly;
- Replace as necessary, at least 1pc per daily.

## Warning

- Single Use only;
- Do not use if package is damaged;
- Keep out of the reach of children.

