



Isolation Gown

Disposable·Waterproof

Isolation Gown provide impervious protection and a comfortable fit. Isolation Gown is manufactured from non-woven material and available in yellow/blue color. Isolation Gowns have universal, full-length size and feature comfortable knitted cuffs.

Applicable for

Economical protective gown for personal and visitors

Specifications

- Item: PP 20gsm+PE 15gsm isolation gown
- Material: PP +PE material,waterproof
- Fabric: PP 20gsm/ m² +PE 15gsm/ m²
- Color: Blue, Yellow
- Size: universal,length 120cm,width 140cm.
- Feature: Knitted Cuffs,Tie back enclosure
- Type: Disposable,Non medical,non sterile



Packing Details

- Package: 10pcs/bag,10bags/ctn
- Carton size: 55x28x30cm
- G.W/N.W: 8.5/7.00kgs

Isolation Gowns standard (According to GB/T 38462-2020)

- FDA Registration Number: 3017234461
- Owner/Operator Number: 10073875
- Proprietary Name: Disposable Non-woven Gown
- Product Code: OEA

