



Face Mask

Disposable

Product Details

Product name: 3ply Medical Face Mask Type IIR
 Product type: Single use
 Size: 17.5cmx 9.5cm
 Layer: 25gsmPP+25gsm meltblown filter+25gsmPP
 Colour: blue or white
 Country of origin: China



Mask performance (according to EN14683:2019+AC:2019 standard)

Bacterial Filtration Efficiency (BFE) :≥ 98%

Instructions



Flatten the mask, and simply bend the nose bar to fit the face shape.



Pull elastic earloops over ears, and extend mask to fully cover mouth, nose and chin.



Mold the nose bar to ensure a snug fit.

Features & Benefits

COMFORTABLE EARLOOP

Comfortable elastic earloop, extra-soft earloops eliminate pressure to ears. You feel comfortable with periods of time usage.

HIDDEN NOSE STRAPS

Inset with bendable nose strips. Freely adjust shape according to face shape.

Storage Conditions

Storage temperature: -20-38°C
 Storage humidity: <80%
 Protected from light, in a dry indoor environment, valid for 5 years.

Logistic Information

Outer case specifications

Size	Gross weight	QTY
52x40x30cm	8.00kg	2000pcs

Inner case specifications

Size	Gross weight	QTY
19x10x7cm	0.2kg	50pcs

