



# Surgical Hand Brush with Chlorhexidine Disposable

Surgical scrub brushes are used to significantly reduce the number of micro-organisms on the hands and forearms before performing surgeries. The sponge on the back of the brush is used for dipping in and sucking disinfecting fluid for cleaning parts rubbed. The nail cleaner feature makes the removal of micro-organisms under the nail quick and easy.

## Specifications

- Item: Surgical Hand Brush with Chlorhexidine 4%
- Spec.: Sponge Chlorhexidine 20ml
- Material: Foam sponge with plastic holder  
Soft brush body made of LDPE
- Content: Foam, nail cleaner, surgical brush
- Function: Disposable
- Sterile: EO gas

## Intended use

Used to Hospital and clinic nursing.

## Packing details

- Package: 250pcs/ctn
- Carton size: 54\*30\*33cm

