



Transport Tube with Amies

Disposable

Amies transport medium is a widely used and effective semisolid medium for the transportation of swab specimens to the microbiology laboratory. Placing swabs in a moist container or transport medium prevents drying and death of bacteria.

Amies is a modification of Stuart's medium in which glycerophosphate is replaced by an inorganic phosphate buffer.

Intended use

Can be used to collect specimens

Specifications

- Item: Transport Tube with Amies
- Tip Material: viscose, rayon, nylon
- Medium Amies
- Length: 12*160 MM
- Stick: PS snappable stick
- Sterile: EO or Gamma



Packing details

- Package: 1pcs*100bag*10bag/ctn
- Carton size: 62*32*45cm
- G.W/N.W: 14.5/13.5kgs

Transport Tube Safety Standard (According to EN 149 -2001+A1-2009)

