



Dental Micro Applicators

Disposable

Disposable Dental Micro Applicators are high-quality micro applicators. Bend into to any angle, especially in areas of limited access. The material to be applied is held in suspension between the fibers, creating a non-drop sphere that retains the material until it is applied.

Applicable for

Used for the application of any kind of liquids, pastes, or gels used in dentistry.

Specifications

·Item:	Dental Micro Applicators
·Material:	PP + Nylon
·Diameter:	2.5mm,2.0mm,1.5mm,1.2mm
·Color:	Blue,Green,Purple,White
·Characteristic:	Bends into any angle.

Features




1. Convenient handle
2. Disposable,non-linting, non-absorbent
3. Colored for easy identification
4. Available Sizes: Regular, Fine, Ultarfine, and Cylinder.





Details & Packing

Blue Green Purple White
Φ Φ Φ Φ
2.5mm 2.0mm 1.5mm 1.2mm



Packing Details

Package: 100pcs/barrel,4barrels/box,40boxes/ctn
Carton size: 51.5*43*25cm
G.W./N.W.: 14/13kgs/ctn

