



Disposable Coverall

Disposable coverall provide impervious protection and a comfortable fit. Disposable coverall is manufactured from non-woven material and available in blue/white color. Disposable coverall have S,M,L,XL,XXL,XXXL, full-length size and feature comfortable elastic cuffs.

Intended use

Economical protective gown for personnal and visitors

Specifications

- Item: Disposable Coverall
- Material: PP
- Fabric: PP 35-60gsm/ m²
- Color: Blue,White
- Size: S,M,L,XL,XXL,XXXL
- Feature: Elastic cuffs.,Tie back enclosure
- Type: Disposable,non sterile



Packing details(35gsm)

- Package: 1pcs/bag,50pcs/ctn
- Carton size: 52x26x45cm
- G.W/N.W: 6.7/6.2kgs

Disposable Coverall standard (According to GB15979-2002)

- FDA Registration Number: 3017234461
- Owner/Operator Number: 10073875
- Proprietary Name: Disposable Non-woven Gown
- Product Code: OEA

