



IV Administration Set

Disposable

IV Administration Sets (iv sets) provide healthcare providers with a fast, effective way to deliver fluids and medication to patients. Infusion IV therapy can be implemented in a variety of situations, such as Dehydration, Electrolyte imbalances, and Specialized medication delivery.

Intended use

Can be used for intravenous administration.

Specifications

- Sterile
- Latex Free
- Without Needle
- Non-Vented Spike
- 10 Drops ml with Guard/Drip Chamber 60mm
- With Y-injection Site
- Available with Luer Lock and Luer Slip Connection Tips
- 3.0mm x 4.1mm PVC Tubing-Length 135+15cm



Packing details

- Package: PEbag,25pcs/bag,500pcs/ctn
- Carton size: 59*38*32.5cm
- G.W/N.W: 12/11kgs

