



## Isolation Gown

## Disposable·PP

Isolation Gown provide impervious protection and a comfortable fit. Isolation Gown is manufactured from non-woven material and available in blue color. Isolation Gowns have universal, full-length size and feature comfortable elastic cuffs.

## Applicable for

Economical protective gown for personal and visitors

## Specifications

- Item: PP 24gsm isolation gown
- Material: PP material
- Fabric: 24gsm/m<sup>2</sup>
- Color: Blue
- Size: universal,length 120cm,width 140cm.
- Feature: Elastic Cuffs,Tie back enclosure
- Type: Disposable,Non medical,non sterile



## Packing Details

- Package: 10pcs/bag,10bags/ctn
- Carton size: 52x26x30cm
- G.W/N.W: 6.5/5.5kgs

## Isolation Gowns standard (According to GB/T 38462-2020)

- FDA Registration Number: 3017234461
- Owner/Operator Number: 10073875
- Proprietary Name: Disposable Non-woven Gown
- Product Code: OEA

