

2way Catheter



3way Catheter


Silicone Foley Catheter

Disposable

The Silicone Foley Catheter is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bedridden.

Specifications

- Item: Foley Catheter
- Material: Silicone
- Length: 260-400mm
- Size: Fr6 - Fr26
- Type: 2way/3way
- Sterile: EO



Features

1. Made from 100% medical grade silicone;
2. Good for long term placement;
3. X-ray detective line through the catheter;
4. Color-coded for visualization of the size;
5. Length: 310mm (pediatric); 400mm (standard).





Specification Details

Model	Size(Fr/Ch)	Balloon(CC)	Color	Length(mm)
2-way pediatric	6	3	Light red	270
	8	5	Black	270
	10	5	Grey	270
2-way Standard	12	15	White	400
	14	15	Green	400
	16	15	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
	26	30	Pink	400
	2-way female	12	15	White
14		15	Green	260
16		15	Orange	260
18		30	Red	260
20		30	Yellow	260
22		30	Violet	260
3-way standard	16	30	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
	26	30	Pink	400

Packing Details

Package: 10pcs*20boxes/ctn

Carton size: 52*34*24.5cm

G.W/N.W: 7/6kgs

