



## Specimen Collection

## Urine Plastic Container

Urine Plastic Containers are individually wrapped and made from medical grade polypropylene material, Each urine cup has a screw-on leak-proof lid that can be sealed tightly for safe, mess-free handling. Resist high temperatures up to 121°C, autoclavable. Various shapes, volume, and color designs for different specimen collection and test requests.

### Intended use

Used for urine collection and testing in laboratories, clinics, and hospitals.

### Features

- 1.Clear molded graduation for reading and large frosted area for marking and writing.
- 2.Good sealing ensures quality control before the inspection.
- 3.Available with customized bar code.
- 4.Available in individual or bulk packs.

### Specifications

· Item:	Urine Plastic Container
· Material:	PP
· Cap Structure:	Screw Cap
· Volume:	120ml
· Sterile:	E.O
· Packing:	400pcs/ctn
· G.W/N.W:	8/6.8kgs

