



Bactericidal Wipes

70% isopropyl alcohol

Bactericidal Wipes is produced by non-woven and 70% isopropyl alcohol which has slightly irritating smell. The product is packaged with the Plastic bucket. One wipes saturated with 70% isopropyl alcohol. Single use only which is safe and continent. The product will not cause adverse stimulation to human skin and mucosa, will not cause allergic reactions and other damage, and will not precipitation toxic substances either.

Intended Use

Used to disinfect hard surface of non-invasive medical device.

Specifications

- Item: Bactericidal Wipes
- Material: Non-woven 35g/ m²
- Content: 70% Isopropyl Alcohol
- Function: Disinfection
- Size: 14x20cm
- Packing method: Plastic bucket





Directions

- One pad saturated with 70% Isopropyl Alcohol.
- Rub skin briskly in a circular motion from injection and injury site outward.
- For External single use only.
- Discard after single use.

Caution

- Not used to disinfect invasive medical device.
- Do not use when the package is unpack or damaged.
- In the event of skin irritation or rash,discontinuance
- Avoid contact with eyes,keep away from children.

Storage and transportation

- The device shall be stored at a relative humidity of less than 80%, the storage environment shall be kept non-corrosive gases, dry, avoiding sunlight, well-ventilated conditions, and clean.
- The device shall be taken care during transportation and handling.

Packing Details

Package:	100pcs/bucket, 20bucket/ctn
Carton size:	51x41x19cm
G.W/N.W:	8.5kgs

