



## Dental Micro Applicators

## Disposable

Disposable Dental Micro Applicators are high-quality micro applicators. Bend into to any angle, especially in areas of limited access. The material to be applied is held in suspension between the fibers, creating a non-drop sphere that retains the material until it is applied.

### Applicable for

Used for the application of any kind of liquids, pastes, or gels used in dentistry.

### Specifications

·Item:	Dental Micro Applicators
·Material:	PP + Nylon
·Diameter:	2.5mm,2.0mm,1.5mm,1.2m
·Color:	m
·Characteristic:	Blue,Green,Purple,White Bends into any angle.

### Features

1. Convenient handle
2. Disposable,non-linting, non-absorbent
3. Colored for easy identification
4. Available Sizes: Regular, Fine, Ultarfine, and Cylinder.





## Details & Packing




Blue      Green      Purple      White

Φ  
2.5mm

Φ  
2.0mm

Φ  
1.5mm

Φ  
1.2mm



## Packing Details

Package: 100pcs/barrel, 4barrels/box, 40boxes/ctn

Carton size: 51.5\*43\*25cm

G.W./N.W.: 14/13kgs/ctn

