



Disposable Coverall

Disposable coverall provide impervious protection and a comfortable fit. Disposable coverall is manufactured from non-woven material and available in blue/white color. Disposable coverall have S,M,L,XL,XXL,XXXL, full-length size and feature comfortable elastic cuffs.

Intended Use

Economical protective gown for personnal and visitors.

Specifications

·Item:	Disposable Coverall
·Material:	PP
·Fabric:	PP 35-60gsm/ m ²
·Color:	Blue,White
·Size:	S,M,L,XL,XXL,XXXL
·Feature:	Elastic cuffs.,Tie back enclosure
·Type:	Disposable,non sterile



Packing Details

Package:	1pcs/bag,50pcs/ctn
Carton size:	52x26x45cm
G.W/N.W:	6.7/6.2kgs

Disposable Coverall standard (According to GB15979-2002)

FDA Registration Number:	3017234461
Owner/Operator Number:	10073875
Proprietary Name:	Disposable Non-woven Gown

