



Face Mask

Disposable

Product Details

Product name: Child Face Mask Earloop 3ply  
Product type: Non- Medical Grade  
Size: 3.5”X5.7”(9cmX14.5CM)  
Inner & Outer layer: 25g Non-woven fabric  
Filter: 25g Meltblown fabric  
Colour: Blue  
Country of origin: China  
Efficiency (BFE) : ≥ 99%

Mask performance (according to GB/T32610-2016 standard)  
Bacterial filtration



Instructions



Flatten the mask,and simply bend the nose bar to fit the face shape.



Pull elastic earloops over ears,and extend mask to fully cover mouth, nose and chin.



Mold the nose bar to ensure a snug fit.

Features&Benefits

COMFORTABLE EARLOOP

Comfortable elastic earloop, extra-soft earloops eliminate pressure to ears. You feel comfortable with periods of time usage.

HIDDEN NOSE STRAPS

Inset with bendable nose strips. Freely adjust shape according to face shape.

SHIELD YOURSELF

Perfect when traveling through congested airports, bus terminals, malls, and busy city streets.

Storage Conditions

Storage temperature:-20-38℃  
Storage humidity:<80%  
Protected from light,in a dry indoor environment,valid for 5years.

Logistic Information



Specifications  
outer case

Size	Gross weight	QTY
52x32x30cm	7.5kg	2000pcs(50pcsX40boxes)



Specifications  
inner case

Size	Gross weight	QTY
15x10x7cm	0.18kg	50pcs

