



Face Mask

Disposable

Product Details

Product name:	3-ply Earloop Non-Woven Face Mask
Product type:	Single use, non-sterile
Size:	17.5cmx9.5cm
Inner & Outer layer:	Non-woven fabric
Filter:	Meltblown fabric
Colour:	Blue/Pink/Yellow/Green/Purple/Black
Country of origin:	China
Efficiency (BFE) :	≥ 95%



Mask performance (according to GB/T32610-2016 standard)
Bacterial filtration

Instructions



Flatten the mask,and simply bend the nose bar to fit the face shape.



Pull elastic earloops over ears,and extend mask to fully cover mouth, nose and chin.



Mold the nose bar to ensure a snug fit.

Features&Benefits

COMFORTABLE EARLOOP

Comfortable elastic earloop, extra-soft earloops eliminate pressure to ears. You feel comfortable with periods of time usage.

HIDDEN NOSE STRAPS

Inset with bendable nose strips. Freely adjust shape according to face shape.

SHIELD YOURSELF

Perfect when traveling through congested airports, bus terminals, malls, and busy city streets.

Storage Conditions

Storage temperature:-20-38℃
Storage humidity:<80%
Protected from light,in a dry indoor environment,valid for 5years.

Logistic Information



Specifications
outer case

Size	Gross weight	QTY
52x38x30cm	8.00kg	2000pcs



Specifications
inner case

Size	Gross weight	QTY
18x10x7cm	0.2kg	50pcs

