



IV Administration Set

Disposable

IV Administration Sets (iv sets) provide healthcare providers with a fast, effective way to deliver fluids and medication to patients. Infusion IV therapy can be implemented in a variety of situations, such as Dehydration, Electrolyte imbalances, and Specialized medication delivery.

Intended Use

Can be used for intravenous administration.

Packing Details

Package:	PEbag, 25pcs/bag, 500pcs/ctn
Carton size:	59*38*32.5cm
G.W/N.W:	12/11kgs



Specifications

·Item	IV Administration Set with Y
·Material	Medical Grade PVC/PP/ABS Regulator/Latex Free
·Spike	Non-Vented Spike
·Needle	Without Needle
Type	10 Drops ml with Guard/Drip Chamber 60mm
·Injection Port	With Y-injection Site
·Connector	Available with Luer Lock and Luer Slip Connection Tips
·Tube Length	3.0mm x 4.1mm PVC Tubing-Length 135+15cm
·Sterilization	EO Gas

