



Infusion Set

Disposable



The disposable infusion set for adult use is a sterile, single-use medical device, applied in venous infusion in the clinic. Equipped with Spike, Air Filter, Flexible tubing, Y-injection set, Roller Clamp, Fluid Filter, needle, etc.

Intended Use

Can be used for intravenous administration.

Specifications

·Item	IV Administration Set with Y
·Material	Infusion tube: PVC Vented Spike: ABS Drip Chamber: PVC Flow regulator: ABS or PE Injection port: Natural latex or latex-free
·Options available	With or without Air vented spike; With or without Hypodermic Needle (19G-27G); With or without Fluid filter; With or without "Y" injection port; Luer lock or Luer slip connector; Or other devices as your request.
·Tube length	150cm or customized
·Drop Type	20 Drops/ML(Pediatric), 60 Drops/ML(Adult)
·Sterilization	Sterile by EO gas , non-toxic , non-pyrogen
·Shelf Life	3 years





Features

1. Equipped with Spike, Air Filter, Flexible tubing, Y-injection set, Roller Clamp, Fluid Filter and .Hypodermic needle, Drip chamber, Two-way connector(luer./luer lock conncetor, etc.
2. Can be assembled with infusion bags and bottles;
3. Sterilized by E.O gas, non-toxic and pyrogen-free;
4. Highly sharp needle points, causing light puncture force and bringing less pain to patient;
5. Flexible tubing with standard length of 150cm; 1.8m, 2.0m or customized length.
6. With flow regulator to control the rate of medicine liquid.

Caution

1. Not used to disinfect invasive medical device.
2. Do not use when the package is unpack or damaged.
3. In the event of skin irritation or rash,discontinuance
4. Avoid contact with eyes,keep away from children.

Storage and Transportation

1. The device shall be stored at a relative humidity of less than 80%, and the storage environment shall be kept non-corrosive gases, dry, avoiding sunlight, well-ventilated conditions, and clean.
2. The device shall be taken care of during transportation and handling.

Packing Details

Package:	PEbag,25pcs/bag,500pcs/ctn
Carton size:	59*38*32.5cm
G.W/N.W:	12/11kgs

