



## Insulin Syringe

## Disposable

Syringes for insulin users are designed for standard U-100 insulin. The dilution of insulin is such that 1 ml of insulin fluid has 100 standards "units" of insulin. Since insulin vials are typically 10 ml, each vial has 1000 units.

### Directions

Tear the single bag of the syringe, remove the syringe with needle, remove the syringe needle protection sleeve, pull the plunger back and forth slide, tighten the injection needle, and then into the liquid, needle up, slowly push the plunger to exclude air, subcutaneous or Intramuscular injection or blood.

### Specifications

·Item:	Insulin syringe
·Material:	Barrel Plunger: Medical Grade PP Needle: Stainless steel Piston: Natural Latex,Synthetic Rubber
·Volume:	0.5ml,1ml
·Needle Gauge:	with integrated needle,27G-30G
·Feature:	Disposable,Medical use



### Packing Details

Package:	100pcs/box, 50boxes/ctn
Carton size:	36x20x17cm
G.W/N.W:	8.5/8kg

