



Isolation Gown

Disposable

Isolation Gown provide impervious protection and a comfortable fit. Isolation Gown is manufactured from non-woven material and available in blue color. Isolation Gowns have universal, full-length size and feature comfortable elastic cuffs.

Intended Use

Economical protective gown for personnel and visitors.

Specifications

- Item: PP 25gsm isolation gown
- Material: PP material
- Fabric: 25gsm/ m²
- Color: Blue
- Size: universal,length 115cm,width 137cm.
- Feature: Elastic Cuffs,Tie back enclosure
- Type: Disposable,Non medical,non sterile



Packing Details

- Package: 10pcs/bag,10bags/ctn
- Carton size: 52x26x30cm
- G.W/N.W: 6.5/5.5kgs

Isolation Gowns standard (According to GB/T 38462-2020)

- FDA Registration Number: 3017234461
- Owner/Operator Number: 10073875
- Proprietary Name: Disposable Non-woven Gown
- Product Code: OEA

