



Isolation Gown

Disposable

Isolation Gown provide impervious protection and a comfortable fit. Isolation Gown is manufactured from SMS material and available in blue color. Isolation Gowns have universal, full-length size and feature comfortable elastic cuffs.

Intended Use

Economical protective gown for personal and visitors.

Specifications

- Item: SMS 42GSM Isolation Gown
- Material: SMS
- Fabric: 42gsm/ m²
- Color: Blue
- Size: universal,length 130cm,width 150cm.
- Feature: Elastic Cuffs,Tie back enclosure
- Type: Disposable,non sterile



Packing Details

- Package: 10pcs/bag,10bags/ctn
- Carton size: 52×38×30cm
- G.W/N.W: 12/11kgs



Isolation Gowns standard (According to GB/T 38462-2020)

- FDA Registration Number: 3017234461
- Owner/Operator Number: 10073875
- Proprietary Name: Disposable Non-woven Gown
- Product Code: OEA

