



2way Foley Catheter



3way Foley Catheter

Latex Foley Catheter

Disposable


Sterile Medical Latex Foley Catheter is made from natural latex, consists of a shaft, drainage funnel, inflation funnel, irrigation funnel (if present), balloon, and valve. This product is for urine drainage use during surgery or post-surgery.

Intended Use

Include providing relief when there is urinary retention, monitoring urine output for critically ill persons, managing urination during surgery, and providing end-of-life care.

Specifications

·Item:	Foley Catheter
·Material:	Latex
·Length:	260-400mm
·Size:	Fr6 - Fr26
·Type:	2way/3way
·Sterile:	EO



Features

1. Made from natural latex;
2. Good biocompatibility;
3. Silicone coated surface reducing allergic reaction;
4. Smooth tapered tip facilitating introduction;
5. Color-coded for visualization of the size;
6. Extra lumen to irrigation and drug delivery;





Specifications Details

Item No	Item	Size	Model	Length	Packing
E20043	Latex Foley Catheter	FR6-Fr26	2-way pediatric	260-400mm	Individual package: Blister/PE packing 10pcs/box,50boxes/ctn
E20044	Silicone Foley Catheter		2-way Standard 2-way female 3-way standard		

Model	Size(Fr/Ch)	Balloon(CC)	Color	Length(mm)
2-way pediatric	6	3	Light red	270
	8	5	Black	270
	10	5	Grey	270
2-way Standard	12	15	White	400
	14	15	Green	400
	16	15	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
	26	30	Pink	400
2-way female	12	15	White	260
	14	15	Green	260
	16	15	Orange	260
	18	30	Red	260
	20	30	Yellow	260
	22	30	Violet	260
3-way standard	16	30	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
	26	30	Pink	400

Packing Details

Package: 10pcs*50boxes/ctn
 Carton size: 52*42*47.5cm
 G.W/N.W: 13/12kgs

