



## PBT Elastic Bandage

## Disposable

PBT elastic bandage is a lightweight bandage with a high level of stretch. The bandages are designed to conform to the body, in turn allowing for freedom of movement. The bandages are lightweight, breathable and fray resistant used together with cast padding.

### Intended Use

Used to provide binding force on wound dressings or limbs to play a role of bandaging and fixation.

### Features

- Lightweight, breathable, and fray resistant;
- Excellent elastic characteristic, soft texture;
- Durable elasticity and stable externals.

### Specifications

·Item:	PBT elastic bandage
·Material:	Composed of Polyester, Spandex, Draw textured yarn
·Fabric gram:	80g/ m <sup>2</sup>
·Size:	15cm*50m
·Color:	White
·Function:	Healthcare
·Glue Type:	Medical Glue or no Glue

