



## Identification Bracelet

## Disposable


Identification bracelet is also called the identification wristband, which is a band worn on the wrist as a mark of personal identification. Mainly used in the medical field, the patient's identity information is recorded on it, which facilitates the hospital's identification and management of the patient, reduces the error rate of medical care, and reduces the incidence of medical accidents.

### Intended Use

For correct Patient Identification

### Specifications

·Item	Identification Bracelet
·Material	Medical Grade PVC
·Color	Pink,blue,white,yellow,red,gray,green
·Type	with card insert write on
·Size	Adult:242*25mm Child:173*25mm Infant:154*25mm
·Closure	Plastic snap for one-time use only
·OEM	Available
·Sample	Available



E20069 Identification Bracelet with card insert




E20070 Identification Bracelet write on





## Specifications Details



Size	Length	B	B1	B2	S	φC	Front	φa φb
Adult	242±2	25±0.2	20±0.2	12.5±0.2	14	φ3.5±0.1	φa=4.0±0.1	
Child	173±2	25±0.2	18±0.2	12.5±0.2	13 14	φ3.5±0.1	φa=4.0±0.1 φb=2.5±0.1	
Infant	154	25±0.2	18±0.2	12.5±0.2	10	φ3.5±0.1	φa=4.0±0.1 φb=2.5±0.1	
Thickness: 0.36-0.38±0.01								


## Packing Details

### Adult Identification Bracelet

Package: 100pcs/box, 50boxes/ctn

Carton size: 56x42x26.5cm

G.W/N.W: 16/14kg




273\*80\*50mm

### Child/Infant Identification Bracelet

Package: 100pcs/box, 50boxes/ctn

Carton size: 42x36x26.5cm

G.W/N.W: 10/9kg



205\*68\*50mm

