



Sample Cup

Disposable

Sample Cup, manufactured from high quality polymer PS material, sample cup for Hitachi Biochemical apparatus.

Features

1. High precision mould makes the cuvettes smooth, bright and clean;
2. Glass-like clarity providing excellent light transmission for testing;
3. Accurate size and structure makes the sample cup adapted to Hitachi biochemical analyzer well.

Specifications

·Item	Sample Cup
·Material	Polystyrene
·Size	•16*38mm

Packing Details

Package:	500pcs/bag,11bags/ctn
Carton size:	56*35*31cm
G.W/N.W:	11.3/10.0kg

