



2way Foley Catheter



3way Foley Catheter

Silicone Foley Catheter

Disposable

The Silicone Foley Catheter is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bedridden.

Intended Use

Include providing relief when there is urinary retention, monitoring urine output for critically ill persons, managing urination during surgery, and providing end-of-life care.

Specifications

| | |
|------------|----------------|
| ·Item: | Foley Catheter |
| ·Material: | Silicone |
| ·Length: | 260-400mm |
| ·Size: | Fr6 - Fr26 |
| ·Type: | 2way/3way |
| ·Sterile: | EO |

Features

1. Made from 100% medical grade silicone;
2. Good for long term placement;
3. X-ray detective line through the catheter;
4. Color-coded for visualization of the size;
5. Length: 310mm (pediatric); 400mm (standard).





Specifications Details

| Item No | Item | Size | Model | Length | Packing |
|---------|-------------------------|----------|--|-----------|--|
| E20043 | Latex Foley Catheter | FR6-Fr26 | 2-way pediatric | 260-400mm | Individual package: Blister/PE packing 10pcs/box,50boxes/ctn |
| E20044 | Silicone Foley Catheter | | 2-way Standard 2-way female 3-way standard | | |

| Model | Size(Fr/Ch) | Balloon(CC) | Color | Length(mm) |
|-----------------|-------------|-------------|-----------|------------|
| 2-way pediatric | 6 | 3 | Light red | 270 |
| | 8 | 5 | Black | 270 |
| | 10 | 5 | Grey | 270 |
| 2-way Standard | 12 | 15 | White | 400 |
| | 14 | 15 | Green | 400 |
| | 16 | 15 | Orange | 400 |
| | 18 | 30 | Red | 400 |
| | 20 | 30 | Yellow | 400 |
| | 22 | 30 | Violet | 400 |
| | 24 | 30 | Blue | 400 |
| | 26 | 30 | Pink | 400 |
| | 12 | 15 | White | 260 |
| | 14 | 15 | Green | 260 |
| 2-way female | 16 | 15 | Orange | 260 |
| | 18 | 30 | Red | 260 |
| | 20 | 30 | Yellow | 260 |
| | 22 | 30 | Violet | 260 |
| | 24 | 30 | Blue | 260 |
| 3-way standard | 16 | 30 | Orange | 400 |
| | 18 | 30 | Red | 400 |
| | 20 | 30 | Yellow | 400 |
| | 22 | 30 | Violet | 400 |
| | 24 | 30 | Blue | 400 |
| | 26 | 30 | Pink | 400 |

Packing Details

Package: 10pcs*50boxes/ctn
 Carton size: 52*42*47.5cm
 G.W/N.W: 13/12kgs

