



Specimen Collection Swab

Disposable

Medical disposable specimen collection swab composed of PP stick and tip. The swab tip material is available in superior absorbent cotton, artificial cotton, and polyester fiber. The tip of the most common applicator is coated with short nylon fibers arranged in a vertical manner.

Intended Use

It is used for sample processing, cell culture, DNA testing, rapid and diagnosis-based testing.

Oral Swab: Oral Sampling

Nasal Swab: Nasal Sampling


Specifications

·Item	Specimen Collection Swab
·Tip Material	Nylon Flocked
·Stick Material	ABS
·Type	Oral Swab,Nasal Swab
·Length	150mm±2mm
·Breakpoint	Nasal Swab: 80mm,Oral Swab:30mm
·Sterilization	EO gas Or without Sterile
·Certification	CE/ISO13485/FDA
·Shelf Life	3 years








Specifications Details



Oral Swab Details



Nasal Swab Details



Packing Details

Package: Individual paper Blister pack, 5000 pcs/ctn

Carton size: 55*28*39cm

G.W/N.W: 9/8kgs

