



Specimen Collection



Sputum Container

Disposable

Sputum containers are individually wrapped and made from medical grade polypropylene(PP) or Polystyrene(PS) material, Each urine cup has a spoon and a screw-on leak-proof lid that can be sealed tightly for safe, mess-free handling. It is used for sputum collection and testing in laboratories, clinics, and hospitals.

Intended Use

Used for sputum collection and testing in laboratories, clinics, and hospitals.

Features

1. Medical grade PP or PS material, ensure reliable quality;
2. Good sealing and watertight, no leakage;
3. Clear molded graduation for reading and large frosted area for marking and writing;
4. Individual poly bag pack, sterile EO gas, no pollution.

Specifications

·Item	Sputun Container
·Material	High Grade PP/PS
·Volume	30ml/40ml/60ml/100ml
·Cap Type	Screw/Hand lift lid
·Cap Color	Orange, Red, Green, White etc
·Cutomize Label	Available
·Package	Bulk or Individual package
·Sterile	With EO sterilization or without Sterile





Specifications Details



Sputum Cup

Urine Cup

Stool Cup

Item No	Type	Material	Cap Structure	Description
E30017	Sputum Container	PP/PS	Screw	30ml/40ml/60ml/100ml
E30018	Urine Container		Screw	30ml/40ml/60ml/100ml/120ml
E30019	Stool Container		Screw	30ml//60ml

Packing Details

Volume	PackingCM	Carton Size(CM)	G.W/N.W(KGS)
30ml	1000pcs/ctn	47*37*59	9/7.8
40ml	800pcs/ctn	47*39*60	8.5/7.5
60ml	600pcs/ctn	56*39*46	7.9/6.6
100ml	400pcs/ctn	60*41*48	8.5/8
120ml	150pcs/ctn	56*40*42	3.5/3

