



## Specimen Collection

## Stool Container

## Disposable

Stool Containers are individually wrapped and made from medical grade polypropylene(PP) material, Each stool container has a screw on leak proof lid that can be sealed tightly for a safe, mess-free handling. It is used for stool collection and testing in laboratory, clinic and hospital.

### Intended Use

Used for stool collection and testing in laboratories, clinics, and hospitals.

### Features

1. Medical grade pp material, ensure reliable quality.
2. Good sealing and watertight, no leakage.
3. Clear molded graduation for read and large frosted area for marking and writing.

### Specifications

·Item	Stool Container
·Material	High Grade PP/PS
·Volume	30ml/40ml/60ml/100ml/120ml
·Cap Type	Screw cover or leakage-proof cap
·Stool	With or without stool
·Color	Cup: Transparent or white Cover: Red or customized
·Package	Bulk or Individual package
·Sterile	With EO sterilization or without Sterile



## Specifications Details



Item No	Type	Material	Cap Structure	Description	Packing
E30017	Sputun Container	PP/PS		30ml	40pcs/bag,25bags/ctn
E30018	Urine Container		Screw	30ml/40ml/60ml/100ml/120ml	16pcs/bag,50bags/ctn
E30019	Stool Container		Screw	30ml//60ml	40pcs/bag,25bags/ctn

## Packing Details

Package: 40pcs/bag,25bags/ctn

Carton size: 47\*37\*59cm

G.W/N.W: 8.9/7.9kgs

