



Transport Tube with Amies

Disposable

Amies transport medium is a widely used and effective semisolid medium for the transportation of swab specimens to the microbiology laboratory. Placing swabs in a moist container or transport medium prevents drying and death of bacteria.

Intended Use

Can be used to collect specimens

Specifications

·Item:	Transport Tube with Amies
·Tip Material:	Viscose, Rayon, Nylon
·Medium	Amies
·Length:	12*160 MM
·Stick:	PS snappable stick
·Sterile:	EO or Gamma

Packing Details

Package:	1pcs*100bag*10bag/ctn
Carton size:	62*32*45cm
G.W/N.W:	14.5/13.5kgs

