



Specimen Collection

Urine Container

Disposable

Urine Plastic Containers are individually wrapped and made from medical grade polypropylene material, Each urine cup has a screw-on leak-proof lid that can be sealed tightly for safe, mess-free handling. Resist high temperatures up to 121°C, autoclavable. Various shapes, volume, and color designs for different specimen collection and test requests.

Intended Use

Used for urine collection and testing in laboratories, clinics, and hospitals.

Features

- 1. Clear molded graduation for reading and large frosted area for marking and writing.
- 2. Good sealing ensures quality control before the inspection.
- 3. Available with customized bar code.
- 4. Available in individual or bulk packs.



Specifications

·Item:	Urine Plastic Container
·Material:	PP
·Cap Structure:	Screw Cap
·Volume:	120ml
·Sterilization:	EO Gas

Packing Details

Package:	1pcs opp bag/box,200pcs/ctn
Carton size:	58*37*35cm
G.W/N.W:	4/3.5kg

