



Disposable Virus Sampling Kit

Disposable

Disposable Virus Sampling Kit: This kit is designed for the collection, transport, and storage of viral samples from oral, nasal, or nasopharyngeal swabs. It includes a sterile swab, a vial of viral transport media (VTM) that inactivates the virus while preserving its RNA, and a biohazard bag.

Intended Use

Ready-to-use system used for the collection and transport of viruses obtained from clinical specimens.

Specifications

·Item:	Disposable Virus Sampling Kit
·Material:	Nylon Flocked tips+ABS stick+Transport Tube
·Liquid Volume:	3ml
·Tube Size:	10ml--16*100mm
·Shape of Tube Bottom:	Skirted V
·Cap Type:	Screw
·Swab:	Nasal Swab,Oral Swab
·Swab Size:	13*150mm
·Shelf Life:	1 Year



Packing Details

Package:	400 kits / ctn
Carton size:	68*31.5*49 cm
G.W/N.W:	8/7 kgs

