



Latex Foley Catheter

Disposable


Sterile Medical Latex Foley Catheter is made from natural latex, consists of a shaft, drainage funnel, inflation funnel, irrigation funnel (if present), balloon, and valve. This product is for urine drainage use during surgery or post-surgery.

Intended Use

Include providing relief when there is urinary retention, monitoring urine output for critically ill persons, managing urination during surgery, and providing end-of-life care.

Specifications

·Item:	Foley Catheter
·Material:	Latex
·Length:	310-400mm
·Size:	Fr6 - Fr24
·Type:	2way/3way
·Sterile:	EO



Features

1. Made from natural latex;
2. Good biocompatibility;
3. Silicone coated surface reducing allergic reaction;
4. Smooth tapered tip facilitating introduction;
5. Color-coded for visualization of the size;
6. Extra lumen to irrigation and drug delivery;





Specifications Details

Item No	Item	Size	Model	Length	Packing
E20043	Latex Foley Catheter	FR6-Fr24	2-way pediatric 2-way Standard 3-way standard	310-400mm	Individual package: Blister/PE packing 10pcs/box,50boxes/ctn

Model	Size(Fr/Ch)	Balloon(CC)	Color	Length(mm)
2-way pediatric	6	3	Light red	310
	8	5	Black	310
	10	5	Grey	310
2-way Standard	12	15	White	400
	14	15	Green	400
	16	15	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
3-way standard	16	30	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400

Packing Details

Package: 10pcs*50boxes/ctn
 Carton size: 52.5x41x43cm
 G.W/N.W: 15/13kgs

