



2way Foley Catheter



3way Foley Catheter

Silicone Foley Catheter

Disposable

The Silicone Foley Catheter is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bedridden.

Intended Use

Include providing relief when there is urinary retention, monitoring urine output for critically ill persons, managing urination during surgery, and providing end-of-life care.

Specifications

- Item: Foley Catheter
- Material: Silicone
- Length: 310-400mm
- Size: Fr6 - Fr24
- Type: 2way/3way
- Sterile: EO



Features

1. Made from 100% medical grade silicone;
2. Good for long term placement;
3. X-ray detective line through the catheter;
4. Color-coded for visualization of the size;
5. Length: 310mm (pediatric); 400mm (standard).





Specifications Details

Item No	Item	Size	Model	Length	Packing
E20044	Silicone Foley Catheter	FR6-Fr24	2-way pediatric 2-way Standard 3-way standard	310-400mm	Individual package: Blister/PE packing 10pcs/box,50boxes/ctn

Model	Size(Fr/Ch)	Balloon(CC)	Color	Length(mm)
2-way pediatric	6	3	Light red	310
	8	5	Black	310
	10	5	Grey	310
2-way Standard	12	15	White	400
	14	15	Green	400
	16	15	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400
3-way standard	16	30	Orange	400
	18	30	Red	400
	20	30	Yellow	400
	22	30	Violet	400
	24	30	Blue	400

Packing Details

Package: 10pcs*50boxes/ctn
 Carton size: 54x52x40cm
 G.W/N.W: 18/16kgs

