



## Specimen Collection

## Urine Plastic Container

## Disposable

Urine Plastic Containers are individually wrapped and made from medical grade polypropylene material, Each urine cup has a screw-on leak-proof lid that can be sealed tightly for safe, mess-free handling. Resist high temperatures up to 121°C, autoclavable. Various shapes, volume, and color designs for different specimen collection and test requests.

### Intended Use

Used for urine collection and testing in laboratories, clinics, and hospitals.

### Features

1. Clear molded graduation for reading and large frosted area for marking and writing.
2. Good sealing ensures quality control before the inspection.
3. Available with customized bar code.
4. Available in individual or bulk packs.



### Specifications

---

·Item:	Urine Plastic Container
·Material:	PP
·Cap Structure:	Screw Cap
·Volume:	120ml
·Sterilization:	EO Gas

---

### Packing Details

---

Package:	1pcs opp bag/box,200pcs/ctn
Carton size:	58*37*35cm
G.W/N.W:	4/3.5kg

---

