



Specimen Collection



Stool Container

Disposable

Stool Containers are individually wrapped and made from medical grade polypropylene(PP) material, Each stool container has a screw on leak proof lid that can be sealed tightly for a safe, mess-free handling. It is used for stool collection and testing in laboratory, clinic and hospital.

Intended Use

Used for stool collection and testing in laboratories, clinics, and hospitals.

Features

1. Medical grade pp material, ensure reliable quality.
2. Good sealing and watertight, no leakage.
3. Clear molded graduation for read and large frosted area for marking and writing.

Specifications

·Item	Stool Container
·Material	High Grade PP/PS
·Volume	30ml/40ml/60ml
·Cap Type	Screw cap
·Stool	With stool
·Color	Cup: Transparent or white Cover: Red or customized
·Package	Bulk or Individual package
·Sterile	With EO sterilization or without Sterile





Specifications Details



Sputum Cup

Urine Cup

Stool Cup

Item No	Type	Material	Cap Structure	Description
E30017	Sputun Container	pp	Hand lift lid/Screw	30ml/40ml/60ml/100ml
E30018	Urine Container		Screw	30ml/40ml/60ml/100ml/120ml
E30019	Stool Container		Screw	30ml/60ml

Packing Details

Volume	PackingCM	Carton Size(CM)	G.W/N.W(KGS)
30ml	40pcs*25bag	47*37*59	9.0/8.5
60ml	1pcs/opp bag, 500pcs/ctn	56*39*46	7.9/6.6

