



Specimen Collection



Urine Container

Disposable

Urine Containers are individually wrapped and made from medical grade polypropylene (PP) or Polystyrene(PS) material, Each urine cup has a screw-on leak-proof lid that can be sealed tightly for safe, mess-free handling. Resist high temperatures up to 121°C, autoclavable. Various shapes, volume, and color designs for different specimen collection and test requests.

Intended Use

Used for urine collection and testing in laboratories, clinics, and hospitals.

Features

- 1. Clear molded graduation for reading and large frosted area for marking and writing.
- 2. Good sealing ensures quality control before the inspection.
- 3. Available with customized bar code.
- 4. Available in individual or bulk packs.

Specifications

| | |
|-----------------|--|
| ·Item | Urine Container |
| ·Material | High Grade PP/PS |
| ·Volume | 30ml/40ml/60ml/100ml/120ml |
| ·Cap Type | Screw cover |
| ·Cap Color | Red,Green,White etc |
| ·Cutomize Label | Available |
| ·Package | Bulk or Individual package |
| ·Sterile | With EO sterilization or without Sterile |





Specifications Details



| Item No | Type | Material | Cap Structure | Description |
|---------|------------------|----------|---------------------|----------------------------|
| E30017 | Sputun Container | PP | Hand lift lid/Screw | 30ml/40ml/60ml/100ml |
| E30018 | Urine Container | | Screw | 30ml/40ml/60ml/100ml/120ml |
| E30019 | Stool Container | | Screw | 30ml/60ml |

Packing Details

| Volume | PackingCM | Carton Size(CM) | G.W/N.W(KGS) |
|--------|------------------|-----------------|--------------|
| 30ml | 1pcs*50bag*20bag | 47*37*59 | 8.5/7.5 |
| 40ml | 1pcs*20bag*40bag | 47*39*60 | 9.0/8.5 |
| 60ml | 1pcs*25bag*20bag | 56*39*46 | 7.9/6.6 |
| 100ml | 1pcs*10bag*30bag | 60*41*48 | 8.5/8 |
| 120ml | 1pcs*10bag*20bag | 58*37*45 | 4.5/4 |

