



Surgical Scalpel Blade with Plastic Handle

Disposable

Surgical Scalpel Blade with Plastic Handle is used for operating at your hospital or surgical center. Surgical Scalpel Blade with Plastic Handle varies by size and shape, while some models are manufactured for specific surgical procedures.

Intended Use

Hospital and clinic, surgery operation.



Specifications

- Item Surgical Scalpel Blade with Plastic Handle
- Material Plastic, Carbon/stainless steel
- Size 6#, 9#, 10#, 11#, 12#, 12B#, 13#, 14#, 15#, 15C#, 16#, 18#, 19#, 20#, 21#, 22#, 23#, 24#, 25#, 26#, 34#, 36#
- Sterilization Gamma Radiation 25kGy
- Certificate CE/ISO13485/FDA
- Shelf Life 5 years
- Feature With blade cover
Non-slip handle
Clearly see the size of each scalp blade
Latex-free







Specifications Details

Item No	Item	Description	Sterile
E20011	Surgical blade	Carbon/stainless steel 6#-36# Gamma Radiation 25kGy	Gamma radiation
E20012	Surgical Scalpel blade with plastic handle	Carbon/stainless steel+Plastic handel 6#-36# Gamma Radiation 25kGy	Gamma radiation
E20012S3	Scalpel handle	Stainless Scalpel handle#3	Gamma radiation
E20012S4	Scalpel handle	Stainless Scalpel handle#4	Gamma radiation

Blade Size

Scalpel Blade Box



Packing Details

Package: 100pcs/box, 50boxes/ctn
 Carton size: 36x20x17cm
 G.W/N.W: 8.5/8kg

Package: 10pcs/box, 50boxes/ctn
 Carton size: 40x32x22cm
 G.W/N.W: 6.5/6kg

