



## Vaginal Speculum

Disposable

The vaginal speculum is a medical device used to dilate and provide access to the vaginal canal for examination and treatment procedures. It facilitates visualization of the cervix and vaginal walls by gently spreading the vaginal walls apart. Vaginal speculums are commonly used during gynecological exams, Pap tests, colposcopy procedures, and certain minor surgical interventions.

### Intended Use

Used for vaginal examination and treatment.

### Specifications

·Item	Vaginal speculum
·Type	Side Screw/Middle screw/normal buckled/ French Type/American Type
·Material	Polystyrene
·Size	Small/Medium/Large
·Package	OPP or PE Individual package
·Sterile	With EO sterilization
·Shelf life	3 years



Type	Material	Size	Package
Vaginal speculum with side screw	Polystyrene	Small/Medium/Large	OPP/PE
Vaginal speculum with middle screw			
Vaginal speculum normal buckled			
Vaginal speculum Fasternet Type(French Type)			
Vaginal speculum Push-Pull Type(America type)			





## Packing Details

Package	Size	Carton Size(CM)	Packing
OPP	S	52x40x37	1pcs/opp bag, 100pcs/carton
	M	52x40x33	1pcs/opp bag, 100pcs/carton
	L	52x40x32	1pcs/opp bag, 100pcs/carton
PE	S	52x39x37	1pcs/pe bag, 100pcs/carton
	M	52x39x33	1pcs/pe bag, 100pcs/carton
	L	52x39x32	1pcs/pe bag, 100pcs/carton



## COA reference

Class	Items of inspection	Technical requirement	Result of inspection	
C	4.2	Size	To be seen standard	1mm
	4.1.6	Tolerance	≤2mm between the head of upper and lower part	Qualified
B	4.1.3	Appearance	Smooth edge from plastic resin-less,	Qualified
			Round and soft shape	
	4.1.5	Strength of attachment	Firm, non-collapsible, non-obstruction	Qualified
	4.1.7	Agility	Flexible, non-obstruction when turning	Qualified
	4.2.1	The resilience distance between the upper and lower part when under 20N pressure	≤30mm	15-16mm
	4.2.2	Strength of parts	Non-broken under 30N pressure when the largest tension	Qualified
A	4.2.1	Sterile	Non-bacteria	Qualified
	4.3.2	Residual of EO gas	≤5μg/g	2.8μg/g

## Caution

1. The speculum is for single use only and must be destroyed after use. second use is prohibited
2. Sterility guaranteed if package unopened. use the speculum immediately after opened.
3. Stored in a cool and dry place.

